

Bioceramic Root Repair

RetroMTA[®]



RetroMTA_Synthesis and hydration behavior of calcium zirconium aluminate (Ca₇ZrAl₆O₁₈) cement

Eun-Hee Kang a, Jun-Sang Yoo b, Bo-Hye Kima, Sung-Woo Choi a, Seong-Hyeon Hong a,*
a Department of Materials Science and Engineering, Research Institute of Advanced Materials, Seoul National University, Seoul 151-744, South Korea
b Department of Oral & Maxillofacial Surgery, School of Dentistry, Seoul National University, Seoul 110-749, South Korea

Perforation Repair

Bactericidal pH 12.5, Calcium hydroxide effect

Sealability Bacterial entombment and fossilization

Biocompatibility New cementum formation

Fast setting 3-min, No-washout & 1-visit treatment


No discoloration Base for esthetics



Package

0.6g Intro Pack 0.3g x 2caps

2.4g Refill Pack 0.3g x 8caps



BioMTA

Seoul National University, Research Park Main Center #323
Gwanak ro1, Gwanak-gu Seoul 151-919, Korea
TEL 02-885-3923 FAX 02-887-3923 <http://www.biomta.com>

Make Teeth Alive